


Fan filter unit (FFU)

Description	Image
<p>Fan filter Unit for Laminar Air flow :</p> <p><u>Size of working FFU size 1225x615x320 mm</u></p> <p><u>HEPA filter size 1220x610x69 mm</u></p> <p>To provide class 100 or ISO class- 5</p> <p>All Stainless Steel 304 Construction Fitted with:</p> <p>HEPA Filter efficiency: 99.999% (No. of 1 HEPA filter)</p> <p>Washable synthetic pre-filter dust holding capacity: 90 %</p> <p>All function control panel: On/Off switches for light & blower.</p> <p>Blower capable to maintain 0.45 m/s face velocity +/- 20%.</p> <p>Variable Speed Control</p> <p>Power: 220V/50Hz</p> <p>GMP compliance model</p> <p>Without Magnehelic Gauge</p> <p>Make: SAF , Hong Kong</p>	

For: Technocratics

Erum Manzoor
0333-0297007