





TECHNOCRATIC SERVICES

ISO 9001:2015 Certified Company

Hollow Cathode Lamp for Atomic Absorption

Description	Name	Photo												
<p>STANDARD HCL FOR ATOMIC ABSORPTION</p> <p>Hollow Cathode Lamp for Atomic Absorption</p> <p>50 mm diameter single-element lamps are designed to be used with our PinAAcle™ and AAnalyst™ Atomic Absorption spectrometer series of instruments.</p> <p>Every genuine PerkinElmer lamp is designed for use with and tested on our spectrometers to assure compatibility and the highest performance.</p> <p>It's ideal for determining most elements by atomic absorption spectroscopy.</p> <p>Specifications</p> <table><tr><td>Element</td><td>Praseodymium</td></tr><tr><td>Model Compatible With</td><td>AAnalyst, PinAAcle</td></tr><tr><td>Number of Elements</td><td>Single</td></tr><tr><td>Outer Diameter</td><td>2 in</td></tr><tr><td>Product Brand Name</td><td>Lumina</td></tr><tr><td>Product Group</td><td>Hollow Cathode Lamp (HCL)</td></tr></table> <p>Make: Lumina, USA</p>	Element	Praseodymium	Model Compatible With	AAnalyst, PinAAcle	Number of Elements	Single	Outer Diameter	2 in	Product Brand Name	Lumina	Product Group	Hollow Cathode Lamp (HCL)	<p>lamp</p> <p>Available lamp</p> <p>Na</p> <p>Ca</p> <p>Hg</p> <p>As</p> <p>Pb</p> <p>Cu</p> <p>Ni</p> <p>Fe</p> <p>etc</p>	 
Element	Praseodymium													
Model Compatible With	AAnalyst, PinAAcle													
Number of Elements	Single													
Outer Diameter	2 in													
Product Brand Name	Lumina													
Product Group	Hollow Cathode Lamp (HCL)													

Head Office Karachi

 +92-21-35001798
 +92-333-0297007
 +92-321-2202241
 +92-300-3972448
 www.technocratics.com.pk
 technoratics@gmail.com
 info@technocratics.com.pk
 Plot No. 763, Future Colony Mor, Main Korangi Industrial Area Road, Karachi-75120

Branch Office Lahore

 0333-0297007, 0335-4944616, 0321-2202241
 Plot # 42, Street # 2, Block-C, Asim Town, Harbanspura Underpass Near Ring Road, Lahore-Pakistan.