





WEIGHING SCALES

Description	Actual picture
<p>Analytical balance Model: TM-EXJ2204H Range: 0.001g to 220g Readibility 0.0001g Pan Ø 118 mm Response Time ≤2s</p> <ul style="list-style-type: none"> • Appearance style hale, fashion and easy, strong motion. • Body made of PC alloy materials, solid and durable. • Scale is made of high quality stainless steel, the modelling simple contracted . • Screen LCD display, backlight is soft and easy to read ; • Equipped with g, ct, oz of 21 kinds of commonly used international quality unit; • Standard RS232 output function, can connect the printer ; • Ac/dc, equipped with battery, rechargeable use ; • Windshield with pure white glass is used to guarantee the stability of the weighing; <p>Rs232 interface for miniprinter Make: Toms, Kong Kong</p>	
<p>Electronic Balance Model: TM-STE3002S Range . 0.2g to 300g Repeatability ±0.01g Pan Ø 118 mm Response Time ≤2s</p> <ul style="list-style-type: none"> • Appearance style hale, fashion and easy, strong motion. • Body made of PC alloy materials, solid and durable. • Scale is made of high quality stainless steel, the modelling simple contracted . • Screen LCD display, backlight is soft and easy to read ; • Equipped with g, ct, oz of 21 kinds of commonly used international quality unit; • Standard RS232 output function, can connect the printer ; <p>Make: Toms, Kong Kong</p>	

Head Office Karachi

+92-21-35001798
 +92-333-0297007
 +92-321-2202241
 +92-300-3972448

www.technocratics.com.pk
 technoratics@gmail.com
 info@technocratics.com.pk

Plot No. 763, Future Colony Mor, Main Korangi Industrial Area Road, Karachi-75120

Branch Office Lahore

0333-0297007, 0335-4944616, 0321-2202241
 Plot # 42, Street # 2, Block-C, Asim Town, Harbanspura Underpass Near Ring Road, Lahore-Pakistan.



TECHNOCRATIC SERVICES

ISO 9001:2015 Certified Company

Electronic Balance

Model: PWT=6KT

Pan size: 225*185mm

Range . 2.00 g to 6.00 kg

Readability 2.0g

- 304 stainless steel attached on pan for contracted and fluent.
- Waterproof scale
- Bigger LCD screen shown elegance and vigorous font.
- With RS232 connected computer, printer, etc , for data transmission
- Portable, AC/DC and lead batteries
- Hole capacity tare and accumulation function, etc.
- Calibration : Internal and External Calibration, Weight Selected

Make: Toms, Kong Kong



Electronic Balance

Model: 3KT

Capacity : 1.0g to 3.00 kg

- Readability 1g
- 304 stainless steel attached on pan for contracted and fluent
- Bigger LCD screen shown elegance and vigorous font.
- Calibration : Internal and External Calibration, Weight Selected
- Portable, AC/DC and lead batteries
- Hole capacity tare and accumulation function, etc.

Make: Toms, Kong Kong



Top Balance from 1GM to 15 KG for Raw Material

Type: Digital

Model: Dove

Display: 20mm RED LED and LCD Display
brightness green backlight

with high

Power: AC adaptor: 9V/500mA: Rechargeable

battery 6v4AH

Pan size: 280x260mm/8.60x10.23 inch Stain

less steel

Net Weight: 4 Kg

Make: Axis Japan



Head Office Karachi



+92-21-35001798



+92-333-0297007



+92-321-2202241



+92-300-3972448

Plot No. 763, Future Colony Mor, Main Korangi Industrial Area Road, Karachi-75120



www.technocratics.com.pk



technoratics@gmail.com



info@technocratics.com.pk

Branch Office Lahore



0333-0297007, 0335-4944616, 0321-2202241



Plot # 42, Street # 2, Block-C, Asim Town, Harbanspura Underpass Near Ring Road, Lahore-Pakistan.