





## STAINLESS STEEL DISTILLATION UNIT FOR LABORATORY

Description	IMAGE
<p><b>Electric Water Single Distilling Apparatus:</b> The pure distilled water is one of the most useful liquids in the Pharmaceutical industries, hospitals, healthcare centers, research institutes and analytic laboratories This still is specially designed for these places to prepare the pure distilled water will be produced through this apparatus as the quality of the raw water supping into the still meets with the specifications of the drinking water of the city.</p> <p>1.electric distilled water machine 2.ISO9001 certificate 3.capacity:5L/h 4.best price</p> <p><b>Technical Specifications :</b> Model Specifications: YAZD-5 Output: 5 liters/hour Power: AC220V/ 4.5KW Package: 540×270×490mm Weight: 9 kg Make: Taiwan</p>	



### Head Office Karachi

+92-21-35001798  
+92-333-0297007  
+92-321-2202241  
+92-300-3972448

www.technocratics.com.pk  
technocratics@gmail.com  
info@technocratics.com.pk

Plot No. 763, Future Colony Mor, Main Korangi Industrial Area Road, Karachi-75120

### Branch Office Lahore

0333-0297007, 0335-4944616, 0321-2202241

Plot # 42, Street # 2, Block-C, Asim Town, Harbanspura Underpass Near Ring Road, Lahore-Pakistan.