




Small Scale double distilled water for Eye drops/Herbal/ Homeopathic drops or small scale R&D (Parental) Production

Description	Picture								
<p>Water Double Distilled Water Production 20 liter per hours Type: Multi-column Water Double Distilled Water made of high quality stainless steel 304 (contact part) and non-contatc part made of stainless steel 304. Twin towers, reasonable structure. Condensing coil structure, the cooling effect is good, large volume of water. Dry automatic cut off the heating power, water to be added back to arrive at the job level and continue heating. Initial steaming and re-distilled water control were set. Copper immersion electric heating pipe, high efficiency, long service life.</p>									
<table border="1"> <thead> <tr> <th>Electric Supply</th> <th>Power</th> <th>Water discharge</th> <th>size)</th> </tr> </thead> <tbody> <tr> <td>380V/50Hz</td> <td>32KW</td> <td>≥20L/h</td> <td>56*87*120</td> </tr> </tbody> </table>		Electric Supply	Power	Water discharge	size)	380V/50Hz	32KW	≥20L/h	56*87*120
Electric Supply		Power	Water discharge	size)					
380V/50Hz		32KW	≥20L/h	56*87*120					
<p>Water quality specification</p>									
<p>Conductivity: 1.3 ms/cm at 25°C</p>									
<p>TOC: <500 ppb</p>									
<p>Total aerobic viable count: NMT 100 CFU/ml</p>									
<p>Total <i>E. coli</i>: Not detectable</p>									
<p>Heterotrophic Count, Not detectable</p>									
<p><i>Pseudomonas aeruginosa</i> Not detectable</p>									
<p><i>Burkholderia cepacia</i> Not detectable</p>									
<p>Sterility 10 CFU/100ml</p>									
<p>Endotoxin content: 0.25 EU/ml</p>									
<p>Make: China</p>									



TECHNOCRATIC SERVICES

ISO 9001:2015 Certified Company



Head Office Karachi



+92-21-35001798



+92-333-0297007



+92-321-2202241

+92-300-3972448



Plot No. 763, Future Colony Mor, Main Korangi Industrial Area Road, Karachi-75120



www.technocratics.com.pk



technocratics@gmail.com

info@technocratics.com.pk

Branch Office Lahore



0333-0297007, 0335-4944616, 0321-2202241



Plot # 42, Street # 2, Block-C, Asim Town, Harbanspura Underpass Near Ring Road, Lahore-Pakistan.