


## LEAK TEST APPARATUS

Description	Image
<p><b>Leak Test Apparatus</b></p> <ul style="list-style-type: none"> <li>The Leak Test Apparatus can be used for quality control purposes to find leaks, holes or other such deficiencies in packaging of tablet blisters packs and strips, bottles, vacuum sealing bags and confectionery products</li> <li>Italian Vacuum pump</li> <li>Oil and maintenance free vacuum pump</li> <li>No return valve system</li> <li>Italian Desiccators PCB based with gas kit</li> <li>Teflon based pressure tolerable internal piping</li> <li>Two vacuum locks LCD display timer</li> <li>Vacuum gauge</li> <li>Stainless steel connectors</li> <li>This is a compact model supporting GMP standards</li> <li><b>Specifications:</b></li> <li>Vacuum level: up to 600 mm Hg</li> <li>Desiccators' sizes: 6, 10, 12 inches</li> <li>Display: Gauge type display</li> <li>Input: Mechanical knob and tactile keyboard</li> <li>Power supply: 220/230 V AC, 50/60 Hz less than 25 VA</li> </ul> <p><b>Make: Dawn ,Taiwan</b></p>	 <p>Model: LT-06</p>

### Head Office Karachi

### Branch Office Lahore