



TECHNOCRATIC SERVICES

ISO 9001:2015 Certified Company

STAINLESS STEEL SCOOP

Description	Image
<p>STAINLESS STEEL SCOOP Material : 304 stainless steel Light duty scoops Good toughness, aging resistance, affordable Complete set : 4 Pcs 200ml = (3" x 4"), 500ml = 3-3/4" x 5-1/2" 750ml = 3-3/4" x 8-1/4" 1000ml = 4-1/3" x 9-1/4"</p>	

Head Office Karachi

- +92-21-35001798
- +92-333-0297007
- +92-321-2202241
- +92-300-3972448
- Plot No. 763, Future Colony Mor, Main Korangi Industrial Area Road, Karachi-75120
- www.technocratics.com.pk
- technoratics@gmail.com
- info@technocratics.com.pk

Branch Office Lahore

- 0333-0297007, 0335-4944616, 0321-2202241
- Plot # 42, Street # 2, Block-C, Asim Town, Harbanspura Underpass Near Ring Road, Lahore-Pakistan.