






TECHNOCRATIC SERVICES

ISO 9001:2015 Certified Company

Friability test apparatus

S. #	Description	Qty																				
01	Friability tester: Model: F2 Compliance: USP/BP Disc: single, transparent Operation: Automatic (control through set time and negative count) RPM: 25 Revolution Make: China																					
02	Friability Tester Model:F8 LCD display Single Drum with quick-change front system Fixed speed (25 rpm) Auto stop on test completion Testing mode: Time Easy discharge of test samples at the end of test into separate collection trays. Manufactured as per USP/BP/EP specifications Power supply: 220V/50Hz/40W Make: China																					
03	FRIABILITY TEST APPARATUS Description: The Portable Tablet Friability Test Instrument Friability Tester consists of a control system, a transmission system, a turntable component, etc., and is controlled by a precision control system consisting of a single-chip microcomputer. The instrument has a reasonable structure, high degree of automation, and control. High schedule, high sensitivity, easy operation and reliable operation. The internal structure is full mechanical gear transmission, no belt, no loss parts, long life. Technical Parameters: <table><tr><td>Turntable cylinder diameter:</td><td>φ286mm</td></tr><tr><td>Turntable cylinder depth:</td><td>39mm</td></tr><tr><td>Sample slip height:</td><td>156mm</td></tr><tr><td>Number of turns of the turntable:</td><td>10-990 lap adjustable (default 100 laps)</td></tr><tr><td>Turntable rotation speed:</td><td>25 rpm</td></tr><tr><td>Turntable speed adjustment accuracy:</td><td>±1 rpm</td></tr><tr><td>Power supply:</td><td>220V/50Hz/40W</td></tr><tr><td>Dimensions:</td><td>370×300×340mm</td></tr><tr><td>Weight:</td><td>15 kg</td></tr><tr><td>Make:</td><td>Lasany , India</td></tr></table>	Turntable cylinder diameter:	φ286mm	Turntable cylinder depth:	39mm	Sample slip height:	156mm	Number of turns of the turntable:	10-990 lap adjustable (default 100 laps)	Turntable rotation speed:	25 rpm	Turntable speed adjustment accuracy:	±1 rpm	Power supply:	220V/50Hz/40W	Dimensions:	370×300×340mm	Weight:	15 kg	Make:	Lasany , India	
Turntable cylinder diameter:	φ286mm																					
Turntable cylinder depth:	39mm																					
Sample slip height:	156mm																					
Number of turns of the turntable:	10-990 lap adjustable (default 100 laps)																					
Turntable rotation speed:	25 rpm																					
Turntable speed adjustment accuracy:	±1 rpm																					
Power supply:	220V/50Hz/40W																					
Dimensions:	370×300×340mm																					
Weight:	15 kg																					
Make:	Lasany , India																					