




# TECHNOCRATIC SERVICES

ISO 9001:2015 Certified Company

## TABLET HARDNESS TESTER

Description	Image
<p><b>Tablet Hardness Tester,</b>  <b>Model YD-1</b>          Tablet hardness testers are instruments for detecting breaking hardness of tablet.</p> <p><b>Features:</b>          The tablet is loaded and pressed artificially.          It has high accuracy press sensor and hardness data is displayed on LED.          It tests tablet hardness continuously.          Test data is displayed and locked automatically, so can test be reset and cycled. Automation: auto-test, auto-diagnose, auto-alarm.</p> <p><b>Technical Specifications:</b>          Hardness range: (2~199.9)N          Hardness accuracy Max: (<math>\pm 0.05\%</math>)          Moving distance of probe: 20mm          Power: 220V/ 50Hz/ 10W or 110V/ 60Hz/ 10W          Dimension: 320*200*100mm  <b>Make:</b> China</p>	

### Head Office Karachi

 +92-21-35001798  
 [www.technocratics.com.pk](http://www.technocratics.com.pk)  
 +92-333-0297007  
 [technoratics@gmail.com](mailto:technoratics@gmail.com)  
 +92-321-2202241  
 [info@technocratics.com.pk](mailto:info@technocratics.com.pk)  
 Plot No. 763, Future Colony Mor, Main Korangi Industrial Area Road, Karachi-75120

### Branch Office Lahore

 0333-0297007, 0335-4944616, 0321-2202241  
 Plot # 42, Street # 2, Block-C, Asim Town, Harbanspura Underpass Near Ring Road, Lahore-Pakistan.