




TECHNOCRATIC SERVICES

ISO 9001:2015 Certified Company

Complete Set-UP of Rotary Evaporator, Pump and Chiller for Herbal Extract & Nutraceutical Extract

Description	Image
<p>Rotary Evaporator (10 Liter collecting Flask) (Herbal & nutraceutical Extraction) Large capacity, large diameter rotary evaporator flask, evaporation area, the vacuum water bath, while rotating edge heating, evaporation of the solution diffusion efficiency can be used for biological, pharmaceutical, chemical, food and other small scale areas, Chinese and production. Products with recycled water multi-use vacuum pumps, diaphragm pumps, low temperature cycle (vacuum) pumps, circulating cooler, constant temperature circulator, low temperature cooling liquid circulating pump and other devices supporting component systems</p> <p>Model RH 1005</p> <p>Rotary Bottle 10L,diameter of flange opening: Φ50mm</p> <p>Collecting bottle 10 L</p> <p>Speed-regulation of main machine Digital display, variable frequency, stepless Speed-regulation</p> <p>Rotary speed of main Machine 5 to 110r/min ,digital display of rotary speed</p> <p>Rotating motor (W) phase AD induction motor,power:40W</p> <p>Condensing tube Vertical type water refrigerating</p> <p>Heating boiler stainless water bath (3KW twin)</p> <p>Rate of Evaporation(H₂O) >3L/h,</p> <p>Temperature control Digital display, temperature control,0~99°C</p> <p>Elevating function Electric elevation</p> <p>Elevating stroke 0~150mm</p> <p>Rated power supply 220V/50Hz</p> <p>Power of power supply 3 KW</p> <p>Make: GlassCo, China</p>	



TECHNOCRATIC SERVICES

ISO 9001:2015 Certified Company

Water Jet Aspirator Pump:

Electrical Requirements: 220V, 50HZ,

Power: **550W**,

Flow: 80L/min, Max vacuum: 0.098Mpa, sucking Speed: 10L/min, Number of suction

Nozzle: 5 units,

Tank Capacity: 34L.

Make: GlassCo, China



Head Office Karachi



+92-21-35001798



+92-333-0297007



+92-321-2202241



+92-300-3972448



Plot No. 763, Future Colony Mor, Main Korangi Industrial Area Road, Karachi-75120



www.technocratics.com.pk



technoratics@gmail.com



info@technocratics.com.pk

Branch Office Lahore



0333-0297007, 0335-4944616, 0321-2202241



Plot # 42, Street # 2, Block-C, Asim Town, Harbanspura Underpass Near Ring Road, Lahore-Pakistan.