




TECHNOCRATIC SERVICES

ISO 9001:2015 Certified Company

Binocular Microscope

Description	Image
<p>Binocular Microscope:</p> <p>Product name: binocular biological microscope</p> <p>Model: XSZ-107BN</p> <p>Total Magnification: 40X~1600X</p> <p>Eyepieces: WF10X(2PCS), WF16X(2PCS)</p> <p>Objectives: 4X, 10X, 40XS, 100XS</p> <p>Eyepieces Tube: Binocular 45°. inclined Tube Rotatable 360°.</p> <p>Illumination: 6V/20W Halogen Lamp</p> <p>Coarse Adjustment Range: 30mm</p> <p>Working Stage: Double Layer Mechanical Stage: 130×140mm</p> <p>Moving Range 75mm×45mm</p> <p>Condenser: NA1.25 Abee Condenser With Changeable Diaphragm, Moving Scale And Fill</p> <p>Size: 32cm*35cm*49cm / 12.6in*13.8in*19.3in</p> <p>Gross Weight: 9 kg/ 20 lb</p> <p>Net Weight: 7.3 kg/ 16.1 lb</p> <p>Make: Hong Kong</p>	

For: Technocratic Service

Erum Manzoor
0333-0297007



Head Office Karachi



+92-21-35001798



+92-333-0297007

+92-321-2202241

+92-300-3972448



Plot No. 763, Future Colony Mor, Main Korangi Industrial Area Road, Karachi-75120



www.technocratics.com.pk



technoratics@gmail.com

info@technocratics.com.pk

Branch Office Lahore



0333-0297007, 0335-4944616, 0321-2202241



Plot # 42, Street # 2, Block-C, Asim Town, Harbanspura Underpass Near Ring Road, Lahore-Pakistan.