







# TECHNOCRATIC SERVICES

ISO 9001:2015 Certified Company

## BIOLOGICAL INDICATOR / CHEMICAL INDICATOR

Sterilization consumable with Japan Pharmacopeial Standard Compliance with ISO 11140/ISO 11138

Description	Image
<p><b>BIOLOGICAL INDICATOR FOR ETHYLENE OXIDE</b></p> <p>This EO Biological indicator are a special type of living microbial products that can be used to confirm the performance of sterilization equipment, the verification of sterilization procedures, and the monitoring of sterilization effects in production processes. It is defined as a defined resistance to a specific sterilization process and is available in the inner packaging for the carrier.</p> <p>Incubation time: 24~48hours</p> <p>Package: 50 pieces/box</p> <p><b>Yellow color means: 123456789 TO</b></p> <p><b>FAILURE (Positive: Purple →Yellow )</b></p> <p><b>Purple color means: 123456789 TO 123456789</b></p> <p><b>SUCCESS (Negative: Purple→ Purple)</b></p> <p><b>Make: Ugaiya, JAPAN</b></p>	 
<p><b>Ethylene Oxide Chemical Indicator Strip:</b></p> <p>This Ethylene Oxide Chemical Indicator Strips are manufactured with special oil/water based heat sensitive ink which leads to chemical reaction accompanied by colour changes under the ethylene oxide sterilization.</p> <p>Size; 95mm×15m</p> <p>Package:200 pieces /Box</p> <p><b>Make: Ugaiya, JAPAN</b></p>	 



# TECHNOCRATIC SERVICES

ISO 9001:2015 Certified Company

## CHEMICAL INDICATOR STRIP FOR STEAM:

The chemical indicator strip is manufactured with special oil/water based heat sensitive ink. The card paper is from special medical grade paper to ensure dryness even in a wet package. The process indicator ink shows clear and accurate color change after a successful sterilization cycle.

Size; 95mm×15mm, PET ,Class4


Packages: 200 pieces/Box

Comply with EN ISO 11140-1

Color changes from yellow to black



**Make: Ugaiya, JAPAN**



UGCI0104, Class4



## BIOLOGICAL INDICATOR FOR STEAM:

These Steam Biological indicators are a special type of living microbial products that can be used to confirm the performance of sterilization equipment, the verification of sterilization procedures, and the monitoring of sterilization effects in production processes. It is defined as a defined resistance to a specific sterilization process and is available in the inner packaging for the carrier.

Incubation time: 24~48hours

Package: 50 pieces/box


**Yellow color means: 123456789 TO**

**FAILURE (Positive: Purple →Yellow )**

**Purple color means: 123456789 TO 123456789**

**SUCCESS (Negative: Purple→ Purple)**

**Make: Ugaiya, JAPAN**



### Head Office Karachi

- +92-21-35001798
- +92-333-0297007
- +92-321-2202241
- +92-300-3972448
- www.technocratics.com.pk
- technocratics@gmail.com
- info@technocratics.com.pk
- Plot No. 763, Future Colony Mor, Main Korangi Industrial Area Road, Karachi-75120

### Branch Office Lahore

- 0333-0297007, 0335-4944616, 0321-2202241
- Plot # 42, Street # 2, Block-C, Asim Town, Harbanspura Underpass Near Ring Road, Lahore-Pakistan.