




TECHNOCRATIC SERVICES

ISO 9001:2015 Certified Company

MINI DRY BLOCK INCUBATOR FOR BIOLOGICAL INDICATOR

Name	Photo
<p>DRY BATH INCUBATOR, DRY BLOCK HEATER, MINI DRY BLOCK INCUBATOR FOR BIOLOGICAL INDICATOR</p> <p>The Block heater (dry bath) is a microprocessor-controlled product . With this technique, more stable and accurate temperature control is achieved. Replacement is convenient with optional metal modules of various specifications configured for different types of test tubes.</p> <p>FEATURES</p> <ol style="list-style-type: none"> 1. LED display 2. Current temperature and count down time display 3. Metal blocks free of sample counter-pollution 4. Easy replacement, cleaning and disinfecting of metal blocks 5. Built-in over-temperature protection 6. Calibrating temperature discrepancy 7. Buzzer alarm when temperature control ends <p>TECHNICAL SPECIFICATIONS:</p> <ol style="list-style-type: none"> 1. Temperature control range: from 5°C above ambient to 100°C 2. Time required for temperature to increase from 20°C to 100°C: ≤15 min 3. Temperature control discrepancy: ≤±0.5°C 4. Display discrepancy: ±0.1°C 5. Homogeneity of block temperature: ≤± 0.5°C 6. Maximum heating power: 150W 7. Maximum temperature: 105°C 8. Time Range: 99 h 59 min 9- Capacity: 21*Φ16mm tubes <p>Make: Hong Kong</p>	

Head Office Karachi



+92-21-35001798



www.technocratics.com.pk



+92-333-0297007



technoratics@gmail.com

+92-321-2202241

info@technocratics.com.pk



Plot No. 763, Future Colony Mor, Main Korangi Industrial Area Road, Karachi-75120

Branch Office Lahore



0333-0297007, 0335-4944616, 0321-2202241



Plot # 42, Street # 2, Block-C, Asim Town, Harbanspura Underpass Near Ring Road, Lahore-Pakistan.