




TECHNOCRATIC SERVICES

ISO 9001:2015 Certified Company

MUFFLE FURNACE WITH BOX TYPE RESISTANCE

Description	Image															
<p><u>MUFFLE FURNACE</u> <u>(box type resistance furnace)</u></p> <p>1 Cold rolling steel spraying exterior; open-side door enables easy on/off. 2. Medium-temperature furnace adopts enclosed fire pot; uniform heating and long service life are provided by spiral electro thermal alloy wire elements that embedded in all four in wall of chamber. 3. High-temperature tubular resistance furnace adopts heat resistant combustion tube and takes elema as heating component to fix on outer sleeve of fire pot. 4. High-temperature box resistance furnace takes elema as heating component fixing in the fire pot to ensure high heat use ratio.</p> <p>Model # SX-5-12 Temperature range: +5 over Ambient -1200 °C Rated power: 5kW</p> <table><tr><td>Packaging dimensions W ×D× H(mm)</td><td>Gross weight (kg)</td><td>Net weight (kg)</td></tr><tr><td>900×590×640</td><td>135</td><td>150</td></tr><tr><td>Interior dimensions W ×D× H(mm)</td><td colspan="2">Exterior dimensions W ×D× H(mm)</td></tr><tr><td>300×200×120</td><td colspan="2">790×500×527</td></tr><tr><td>Power supply</td><td colspan="2">220V</td></tr></table> <p>: Make: TAISITE INSTRUMENT CO.,LTD. China</p>	Packaging dimensions W ×D× H(mm)	Gross weight (kg)	Net weight (kg)	900×590×640	135	150	Interior dimensions W ×D× H(mm)	Exterior dimensions W ×D× H(mm)		300×200×120	790×500×527		Power supply	220V		
Packaging dimensions W ×D× H(mm)	Gross weight (kg)	Net weight (kg)														
900×590×640	135	150														
Interior dimensions W ×D× H(mm)	Exterior dimensions W ×D× H(mm)															
300×200×120	790×500×527															
Power supply	220V															

Head Office Karachi



+92-21-35001798



www.technocratics.com.pk



+92-333-0297007



technoratics@gmail.com
info@technocratics.com.pk



Plot No. 763, Future Colony Mor, Main Korangi
Industrial Area Road, Karachi-75120

Branch Office Lahore



0333-0297007, 0335-4944616, 0321-2202241



Plot # 42, Street # 2, Block-C, Asim Town,
Harbanspura Underpass Near Ring Road,
Lahore-Pakistan.