

Hot plate with magnetic stirrer

<u>Description</u>			Image
<u>Hot plate with stirrer</u>			
Function	Mode	Digital+ Heating + Stirring	
	Max. Volume (ML)	5000	
	Speed (r.p.m)	- 200~2000	
	Working time	Continuous	
	Minimum revolution R.P.M.	200	
	Max. temp(°C)	350	
Structure	Exterior	Iron plate welding, surface spraying	
	Heating plate	Ceramic Coated	
	Driving mode	- Motor	
	Motor power (W)	20	
	Heating power(kW)	0.6	
	Stir bar	PTFE	
Controller	Continuous operation	yes	
	Setting mode	Digital	
	Speed display mode	Digital	
	Timer	no	
	Operation function	no	
	Sensor	PT100	
Specification	Exterior size	320×190×125(W*L*H)mm	
	Packing	390×240×230 W*L*Hmm	
	Plate size mm	190×190	
	Plate load	10KG	
	Power supply50/60 Hz	110V or 220V	
	NW/GW (kg) Make:	5/6 A&E Lab UK Co, Ltd	