



TECHNOCRATIC SERVICES

ISO 9001:2015 Certified Company

MELTING POINT MANUAL

Description	Image
<p><u>Melting Point Tester:</u></p> <p>Model GMP-L1</p> <p>Consist with aluminum cylindrical block for heating, Radial holes admits a thermometer and glass capillary</p> <p>Temp. Controlled by energy regulator</p> <p>Regulator fitted to the unit and with high, low heat switch</p> <p>Temperature range ambient to 350 °C, Resolution 1 /0.1 °C, Freeze the reading, cooling facility to become ready quickly for next test. 220 vac. Max Temp. up to 350°C.</p> <p>Make: India</p>	