



TECHNOCRATIC SERVICES

ISO 9001:2015 Certified Company

Digital LCD Display Microscope:

S. #	Description	Image
01	DIGITAL MICROSCOPE Model: XSZ-100SP Achromatic Objectives: 4X, 10X, 40X (S), 100X (S, Oil) Nosepiece: Quadruple Stage: Double Layer Mechanical Stage 130x125mm Focusing: Coaxial Coarse and Fine Adjustment, Focusing Range 30mm Illumination, condenser NA=1.25:6V/20W Halogen Lamp, Adjustable Brightness Power : AC Input 220V or 110V, DC Output 6V Sensor Resolution: 3.1 Mega Pixel Maximum Resolution: 8 Mega Pixel LCD Screen Size: 5.6 inch LTPS Memory Card: SD/MMC Card USD interface: YES TV output: NTSC/PAL Controller: Remote and USB Shutte Power Supply: DC 12V 1.5A Picture Format: JPEG Video Format: MPEG4 Priter Support: PicBridge and DPDF	

Head Office Karachi



+92-21-35001798



www.technocratics.com.pk



+92-333-0297007



technoratics@gmail.com



+92-321-2202241



info@technocratics.com.pk

Plot No. 763, Future Colony Mor, Main Korangi
Industrial Area Road, Karachi-75120

Branch Office Lahore



0333-0297007, 0335-4944616, 0321-2202241



Plot # 42, Street # 2, Block-C, Asim Town,
Harbanspura Underpass Near Ring Road,
Lahore-Pakistan.