




**TECHNOCRATIC SERVICES**

ISO 9001:2015 Certified Company

## BINOCULAR MICROSCOPES

Description	Image
<p><b>BINOCULAR MICROSCOPES</b></p> <p><b>Model: LB-270</b></p> <p><b>Optic</b> Infinite optical system</p> <p><b>Viewing Head</b> Side and top Binocular viewing head, inclined at 30 degree, 360 degree rotatable, Inter pupillary distance 50-75 mm.</p> <p><b>Eyepiece</b> WF 10X/20 with diopter adjustment WF 16X/13 (optional)</p> <p><b>Objective</b> Infinity Plan Achromatic objectives 4X, 10X, 40X (S), 100X (Oil) (S) (Plan 20X &amp; Plan 60X(S) are optional)</p> <p><b>Nosepiece</b> Backward Quadruple Nosepiece Backward Quintuple Nosepiece (optional)</p> <p><b>Focusing</b> Coaxial coarse and Fine focusing knobs, Travel range: 26mm, scale: 2um</p> <p><b>Condenser</b> Abbe condenser NA 1.25 with Iris Diaphragm</p> <p><b>Filter</b> Blue (Green and yellow filters are optional.)</p> <p><b>Dimension</b> 42cmx28cmx45cm, 8Kg</p> <p><b>Make:</b> Labomed, Inc. USA</p>	



# TECHNOCRATIC SERVICES

ISO 9001:2015 Certified Company

## Head Office Karachi



+92-21-35001798



+92-333-0297007



+92-321-2202241

+92-300-3972448



[www.technocratics.com.pk](http://www.technocratics.com.pk)



[technoratics@gmail.com](mailto:technoratics@gmail.com)

[info@technocratics.com.pk](mailto:info@technocratics.com.pk)

Plot No. 763, Future Colony Mor, Main Korangi  
Industrial Area Road, Karachi-75120

## Branch Office Lahore



0333-0297007, 0335-4944616, 0321-2202241



Plot # 42, Street # 2, Block-C, Asim Town,  
Harbanspura Underpass Near Ring Road,  
Lahore-Pakistan.