




TECHNOCRATIC SERVICES

ISO 9001:2015 Certified Company

MINI DRY BLOCK INCUBATOR FOR BIOLOGICAL INDICATOR

Name	Photo																
<p>MINI DRY BLOCK INCUBATOR FOR BIOLOGICAL INDICATOR: (Portable small incubator)</p> <p>Features:</p> <ul style="list-style-type: none"> - LED display for easy reading - Heating lid - Temperature is up to 100°C - Simultaneous display of set/actual time and temperature - 3D digital wheel adjustment of speed and temperature - High temperature precision is up to $\pm 0.5^{\circ}\text{C}$ - Programmable operation can store up to 5 programs - Low operation voltage (12v) can be used directly with car charger - Various applications with different blocks <p>Technical Specifications:</p> <table> <tr> <td>Temperature Control Range</td><td>RT + 5°C~100°C</td></tr> <tr> <td>Timing Range</td><td>0~99h59min</td></tr> <tr> <td>Temp. Control Accuracy</td><td>$\pm 0.5^{\circ}\text{C}$</td></tr> <tr> <td>Display Accuracy</td><td>$\pm 0.1^{\circ}\text{C}$</td></tr> <tr> <td>Heating Time (25°C to Max. Temp.)</td><td>$\leq 15\text{min}$</td></tr> <tr> <td>1-Heating Lid</td><td>yes</td></tr> <tr> <td>Dimension(L x W x H) mm</td><td>164 x 116 x 98</td></tr> <tr> <td>Net Weight</td><td>0.90kg</td></tr> </table> <p>Make: Mulab, Korea</p>	Temperature Control Range	RT + 5°C~100°C	Timing Range	0~99h59min	Temp. Control Accuracy	$\pm 0.5^{\circ}\text{C}$	Display Accuracy	$\pm 0.1^{\circ}\text{C}$	Heating Time (25°C to Max. Temp.)	$\leq 15\text{min}$	1-Heating Lid	yes	Dimension(L x W x H) mm	164 x 116 x 98	Net Weight	0.90kg	
Temperature Control Range	RT + 5°C~100°C																
Timing Range	0~99h59min																
Temp. Control Accuracy	$\pm 0.5^{\circ}\text{C}$																
Display Accuracy	$\pm 0.1^{\circ}\text{C}$																
Heating Time (25°C to Max. Temp.)	$\leq 15\text{min}$																
1-Heating Lid	yes																
Dimension(L x W x H) mm	164 x 116 x 98																
Net Weight	0.90kg																

Different Standard Power Line

Power Adapter



Head Office Karachi

Branch Office Lahore

+92-21-35001798
 +92-333-0297007
 +92-321-2202241
 +92-300-3972448
 Plot No. 763, Future Colony Mor, Main Korangi Industrial Area Road, Karachi-75120

0333-0297007, 0335-4944616, 0321-2202241
 Plot # 42, Street # 2, Block-C, Asim Town, Harbanspura Underpass Near Ring Road, Lahore-Pakistan.