



TECHNOCRATIC SERVICES

ISO 9001:2015 Certified Company

VORTEX MIXER FOR TESTUBE & BEAKERS

Description	Image
<p>Vortex Mixer for beaker Model: XH-C</p> <p>Operation method: Circling</p> <p>Circling diameter: 3mm / 4mm</p> <p>Rating speed: 2800rpm</p> <p>Working plate size: 55mm</p> <p>Input voltage(frequency): AC 220V±10%(50/60Hz)</p> <p>Permissible relative moisture: 80%</p> <p>Output power: 30W</p> <p>Protective class: IP21</p> <p>Permissible ambient temperature: 5-50°C</p> <p>Net. Weight: 2.5Kg</p>	
<p>Vortex Mixer for test tube Model: XH-D</p> <p>Operation method: Circling</p> <p>Rating speed: 2800rpm</p> <p>External size: 70*120*170mm</p> <p>Working plate: Bowl shape, disc shape</p> <p>Power rating: 50W</p> <p>Working method: Continuous, speed adjustable</p>	