




TECHNOCRATIC SERVICES

ISO 9001:2015 Certified Company

AUTOMATIC WIRE LOOP STERILIZER Bacti-Cinerator Sterilizer

Details	Image														
<p>Bacti-Cinerator Sterilizer sterilizes micro-organisms utilizing infrared heat produced by a ceramic core element. The ceramic element contains no asbestos and ensures maximum sterility without splatter across the work surface area. Complete sterilization occurs within 5-7 seconds at optimum sterilizing temperature of 1500°F (815.6°C).</p> <p>The small footprint makes the Bacti-Cinerator Sterilizer an ideal instrument for anaerobic and aerobic chambers alike. Several user friendly accessories are also available for the Bacti-Cinerator sterilizer.</p> <p>FEATURES:</p> <ul style="list-style-type: none"> The principle of infrared, products are widely used in loop or inoculating needle sterilization, completely replace the alcohol lamp. The temperature of heating hole can reach 800°C, sterilization only takes 5-7 seconds, thorough sterilization. Products can be used in anaerobic chamber. Small volume, light weight, appearance and easy to clean. Heating body can adjust angle, easy operation. Organic substances in the ceramic funnel pipe depths ashing prevent infectious spatter and cross contamination. Precise temperature control technology, avoid the heating components work overload, long life. <table> <tr> <td>Center Max. Temp.</td><td>825°C</td></tr> <tr> <td>Low Level Temp.</td><td>480°C</td></tr> <tr> <td>Max. Sterilize Diameter</td><td>Φ14mm</td></tr> <tr> <td>Heating Element Length</td><td>140mm</td></tr> <tr> <td>Adjustable Angle</td><td>75°~120°</td></tr> <tr> <td>Voltage</td><td>AC 220V , 50Hz</td></tr> <tr> <td>Power</td><td>230W</td></tr> </table> <p>Make: Mulab, Korea</p>	Center Max. Temp.	825°C	Low Level Temp.	480°C	Max. Sterilize Diameter	Φ14mm	Heating Element Length	140mm	Adjustable Angle	75°~120°	Voltage	AC 220V , 50Hz	Power	230W	
Center Max. Temp.	825°C														
Low Level Temp.	480°C														
Max. Sterilize Diameter	Φ14mm														
Heating Element Length	140mm														
Adjustable Angle	75°~120°														
Voltage	AC 220V , 50Hz														
Power	230W														

Head Office Karachi



+92-21-35001798



+92-333-0297007



+92-321-2202241



+92-300-3972448

Plot No. 763, Future Colony Mor, Main Korangi
Industrial Area Road, Karachi-75120



www.technocratics.com.pk



technoratics@gmail.com



info@technocratics.com.pk

Branch Office Lahore



0333-0297007, 0335-4944616, 0321-2202241



Plot # 42, Street # 2, Block-C, Asim Town,
Harbanspura Underpass Near Ring Road,
Lahore-Pakistan.