


## HEATING MANTLE FOR VARIOUS LIQUIDS

Description	Image
<p><b>Heating Mantle:</b> It applies to colleges, petrochemical industries, medicine, and environment protection institutions for liquid heating.</p> <p><b>Characteristics:</b></p> <ol style="list-style-type: none"> <li>1. Cold rolling sheet stretching and spraying exterior.</li> <li>2. Electronic temperature regulation mode with characteristics of stepless speed adjustment, uniform temperature, fast heating-up, and safety and durability.</li> <li>3. Featured by anti-corrosion, age-resistant, durable and safe.</li> <li>4. Beautifully shaped, easy to operate and effective.</li> </ol> <p><b>Specification:</b></p> <p>Model: 98-I-B Capacity (ml): 1Ltr Voltage(V): 220/50Hz Max. Temp. (° C): 450 Power(W): 80 Working time: continuous Exterior Size(mm): Φ 220× 160 Make: China</p>	

### Head Office Karachi



+92-21-35001798



www.technocratics.com.pk



+92-333-0297007



technoratics@gmail.com



+92-321-2202241  
+92-300-3972448  
Plot No. 763, Future Colony Mor, Main Korangi  
Industrial Area Road, Karachi-75120



info@technocratics.com.pk

### Branch Office Lahore



0333-0297007, 0335-4944616, 0321-2202241



Plot # 42, Street # 2, Block-C, Asim Town,  
Harbanspura Underpass Near Ring Road,  
Lahore-Pakistan.