


DRYING OVEN FOR LAB

Description	Image																				
<p>Drying Oven: Capacity: 53Liter</p> <p>Electrical data: 220V 50HZ Temperature range: RT+10-200c Temperature fluctuation: $\pm 1\%$ Power: 800W Chamber dimension: 415X370X345 External dimension (mm): 720X500X510 Locating Range: 1-9999minut Make: China</p>																					
<p>DIGITAL Drying Oven, (220°C) Size: 52Liter Material :Inner stainless steel and outter MS powder coated Several sets of data shows in a LED screen or big LCD screen (when you order, please indicate). The operation interface can be easily understood and handling. - Unique air duct system, no manual adjustment any more. Vapor can be released automatically. - Microcomputer control with timing. Overtemperature deviation protection and digital display, accurate and reliable. - Specular stainless-steel chamber, semicircular arcs at corners for easy cleaning, and the space between the shelves can be set by the customer.</p> <table border="1"> <tr> <td>Work voltage</td><td>AC220V 50HZ</td></tr> <tr> <td>Range of temperature control</td><td>RT+10~280°C</td></tr> <tr> <td>Wave degree of temperature</td><td>+/-1°C</td></tr> <tr> <td>Temperature resolution</td><td>0.1°C</td></tr> <tr> <td>Input power</td><td>850W</td></tr> <tr> <td>Interior dimension W*D*H (mm)</td><td>340*330*320</td></tr> <tr> <td>Overall dimension W*D*H (mm)</td><td>625*540*500</td></tr> <tr> <td>Capacity</td><td>30L</td></tr> <tr> <td>Shelves per chamber (standard equipped)</td><td>2pcs</td></tr> <tr> <td>Timing range</td><td>1~9999min</td></tr> </table> <p>Make: China</p>	Work voltage	AC220V 50HZ	Range of temperature control	RT+10~280°C	Wave degree of temperature	+/-1°C	Temperature resolution	0.1°C	Input power	850W	Interior dimension W*D*H (mm)	340*330*320	Overall dimension W*D*H (mm)	625*540*500	Capacity	30L	Shelves per chamber (standard equipped)	2pcs	Timing range	1~9999min	
Work voltage	AC220V 50HZ																				
Range of temperature control	RT+10~280°C																				
Wave degree of temperature	+/-1°C																				
Temperature resolution	0.1°C																				
Input power	850W																				
Interior dimension W*D*H (mm)	340*330*320																				
Overall dimension W*D*H (mm)	625*540*500																				
Capacity	30L																				
Shelves per chamber (standard equipped)	2pcs																				
Timing range	1~9999min																				

Head Office Karachi

Branch Office Lahore