



TECHNOCRATIC SERVICES

ISO 9001:2015 Certified Company

Standard weights for QC and production 1mg-10kg

Description						Image
STANDARD WEIGHTS SERIES: Material special stainless steel .class F1 Nominal value ; 1mg-10kg No. of Pieces: 29						
Milligram (mg)	Qty	Gram (g)	Qty	Kilo Gram (Kg)	Qty	
1mg	1	1g	1	1Kg	1	
2mg	2	2g	2	2Kg	2	
5mg	1	5g	1	5Kg	1	
10mg	1	10g	1	10Kg	1	
20mg	2	20g	2			
50mg	1	50g	1			
100mg	1	100g	1			
200mg	2	200g	2			
500mg	1	500g	1			
Density 7.96 g/cm Susceptibility.≤0 .01 Magnetization: F1≤ 25μ T Package: Aluminum BOX. Approved by FDA and SGS Standard weight series with calibration certificate. Inlay: From UK. Approved by FDA and SGS Make: Hong Kong						

Please feel free to contact Ms .Erum Manzoor, 021-35001798, 0333-0297007,