



TECHNOCRATIC SERVICES

ISO 9001:2015 Certified Company

Standard weights for QC and production 1mg-200g

Description				Image
<u>STANDARD SERIES:</u> Material special stainless steel class F1 Nominal value ; 1mg-200g Density 7.94-7.96 g/cm Susceptibility. ≤ 0.01 Magnetization: $F1 \leq 25 \mu T$ Package: Aluminum BOX. Approved by FDA and SGS Inlay : From UK. Approved by FDA and SGS No.of pieces: 23				
Miligram (mg)	Qty	Gram (g)	Qty	
1mg	1	1g	1	
2mg	2	2g	2	
5mg	1	5g	1	
10mg	1	10g	1	
20mg	2	20g	2	
50mg	1	50g	1	
100mg	1	100g	1	
200mg	2	200g	2	
500mg	1			

Please feel free to contact Ms .Erum Manzoor, 021-35001798, 0333-0297007,