





TECHNOCRATIC SERVICES

ISO 9001:2015 Certified Company

Overhead Stirrer (SCH)

Description	Actual picture
<p>Overhead Stirrer Model : (SCH)</p> <p>Overhead stirrer is suitable for mixture of high viscosity liquid or solid-liquid.</p> <p>Overhead stirrer is suitable for mixture of high viscosity liquid or solid-liquid. Widely used in chemical synthesis, pharmaceutical, physical and chemical analysis, petrochemical, cosmetics, health care, food, biotechnology and other fields. One year Warranty.</p> <p>Specifications:</p> <p>Max. Volume(L): 20 Speed(r.p.m): 100~2500 Working time: Continuous Minimum revolution (r.p.m): 100 Motor power (W): 60 Setting mode: Digital Speed display mode: Digital Timer: Yes Operation function: Yes Exterior size (W*L*H)(mm): 115×95×240 Packing size (W*L*H)(mm): 400×250×210 Plate size (mm): 200×300 Plate load: 50KG Power supply: 100~240V NW/GW (kg): 6/7</p> <p>Make : A&E Lab Uk, Co. Ltd</p>	 

Please feel free to contact **Ms Erum Manzoor** **021-35001798**, 0333-0297007

Head Office Karachi

+92-21-35001798
 +92-333-0297007
 +92-321-2202241
 +92-300-3972448
www.technocratics.com.pk
technoratics@gmail.com
info@technocratics.com.pk
 Plot No. 763, Future Colony Mor, Main Korangi
Industrial Area Road, Karachi-75120

Branch Office Lahore

+ 92 301 9228811, + 92 334 6600557
 L-32 Main Street, Al-Nafy Homes, Ring Road,
Near WAPDA Office Harbanspura, Lahore