



Autoclave (Bench top) for small operation theater , dental clinic scientific research or microbiology lab

Name	Photo																						
<p>Autoclave (Bench top) /Steam sterilizer</p> <p>Type : horizontal/ front door loading</p> <p>Capacity : 20L(φ250×420 mm)</p> <p>Autoclave is ideal equipment for dental clinics, scientific research or microbiology lab and other organizations to sterilize surgical instruments, fabric glassed, culture media etc</p> <p>Features:</p> <ol style="list-style-type: none"> 1. Indicator light indicates working state. 2. For 4~6 minutes rapidly sterilizing. 3. Sterilizing temperature and time can be preset. 4. Steam-water inner circulation system: no steam discharge, and the environment for sterilizing will be clean and dry. 5. Over temperature & over pressure auto-protection. 6. Safe protection of water lacking. 7. Exhaust the cool air automatically. 8. With three stainless steel sterilizing plates. 9. The chamber of the sterilizer is made of stainless steel. 10. Automatically shut off with beep reminding after sterilization. 	 <p style="text-align: center;">Ready stock</p>																						
<p>Specification</p> <table border="1" data-bbox="113 1299 935 1749"> <tbody> <tr> <td>Sterilizing chamber volume</td> <td>20L(φ250×420 mm)</td> </tr> <tr> <td>Maximum working pressure</td> <td>0.22Mpa</td> </tr> <tr> <td>Maximum working temperature</td> <td>134°C</td> </tr> <tr> <td>Adjustment of temperature</td> <td>105-134°C</td> </tr> <tr> <td>Timer</td> <td>0-60 min</td> </tr> <tr> <td>Chamber temperature equal</td> <td>≤ ± 1°C</td> </tr> <tr> <td>Source power</td> <td>1.5KW/AC 220V 50Hz</td> </tr> <tr> <td>Sterilizing plate</td> <td>340×200×30 mm (3 piece)</td> </tr> <tr> <td>Dimension</td> <td>480×480×384 mm</td> </tr> <tr> <td>Package dimension</td> <td>700×580×500 mm</td> </tr> <tr> <td>G. W/N.W</td> <td>44/40 kg</td> </tr> </tbody> </table> <p>Make: TM. Taiwan</p>	Sterilizing chamber volume	20L(φ250×420 mm)	Maximum working pressure	0.22Mpa	Maximum working temperature	134°C	Adjustment of temperature	105-134°C	Timer	0-60 min	Chamber temperature equal	≤ ± 1°C	Source power	1.5KW/AC 220V 50Hz	Sterilizing plate	340×200×30 mm (3 piece)	Dimension	480×480×384 mm	Package dimension	700×580×500 mm	G. W/N.W	44/40 kg	 
Sterilizing chamber volume	20L(φ250×420 mm)																						
Maximum working pressure	0.22Mpa																						
Maximum working temperature	134°C																						
Adjustment of temperature	105-134°C																						
Timer	0-60 min																						
Chamber temperature equal	≤ ± 1°C																						
Source power	1.5KW/AC 220V 50Hz																						
Sterilizing plate	340×200×30 mm (3 piece)																						
Dimension	480×480×384 mm																						
Package dimension	700×580×500 mm																						
G. W/N.W	44/40 kg																						



TECHNOCRATIC SERVICES

ISO 9001:2015 Certified Company

Product Introduction



Product Details



SUS304 stainless steel heating Pipe--for heating water to produce steam.



Head Office Karachi



+92-21-35001798



+92-333-0297007



+92-321-2202241



+92-300-3972448



Plot No. 763, Future Colony Mor, Main Korangi Industrial Area Road, Karachi-75120



www.technocratics.com.pk



technoratics@gmail.com
info@technocratics.com.pk

Branch Office Lahore



0333-0297007, 0335-4944616, 0321-2202241



Plot # 42, Street # 2, Block-C, Asim Town, Harbanspura Underpass Near Ring Road, Lahore-Pakistan.