






# TECHNOCRATIC SERVICES

ISO 9001:2015 Certified Company

## Low Speed Centrifuge

Description	Actual picture																						
<p><b>Low Speed Centrifuge</b></p> <p><b>Model: TD4Y</b></p> <p>TD4Y low speed centrifuge is portable design and has low noise. It is popular with clinic, research and university.</p> <p><b>Features:</b></p> <ul style="list-style-type: none"> <li>◆ Microcomputer control, brushless motor.</li> <li>◆ It is programmable and has LED display.</li> <li>◆ It has speed/RCF switch easily.</li> <li>◆ There are lid safety device, speed safety device, errors detection.</li> <li>◆ ABS housing make the clean easy</li> </ul> <p>Silastic sealed ring comfort to GMP US FDA certificated</p> <p><b>Technical Specification:</b></p> <table border="1" data-bbox="178 1102 974 1690"> <thead> <tr> <th>Model</th> <th>AE-TD4Y</th> </tr> </thead> <tbody> <tr> <td>Max Speed</td> <td>4000r/min</td> </tr> <tr> <td>Max RCF</td> <td>1980×g</td> </tr> <tr> <td>Max Capacity</td> <td>15mlx6</td> </tr> <tr> <td>Speed Accuracy</td> <td>±30r/min</td> </tr> <tr> <td>Time Setting</td> <td>1min-99min</td> </tr> <tr> <td>Noise</td> <td>&lt;65dB(A)</td> </tr> <tr> <td>Power Supply</td> <td>110V/60Hz, 220V/50Hz</td> </tr> <tr> <td>Dimension(LxWxH)</td> <td>310mm×260mm×220mm</td> </tr> <tr> <td>Packing Size(LxWxH)</td> <td>370mm×320mm×280mm</td> </tr> <tr> <td>Net Weight</td> <td>14kg</td> </tr> </tbody> </table> <p><b>Make :</b> A&amp;E LAB UK Co, Ltd</p>	Model	AE-TD4Y	Max Speed	4000r/min	Max RCF	1980×g	Max Capacity	15mlx6	Speed Accuracy	±30r/min	Time Setting	1min-99min	Noise	<65dB(A)	Power Supply	110V/60Hz, 220V/50Hz	Dimension(LxWxH)	310mm×260mm×220mm	Packing Size(LxWxH)	370mm×320mm×280mm	Net Weight	14kg	  <p>Rotor Capacity: 7ml×12 Max Speed: 4000r/min Max RCF: 1980×g</p>  <p>Rotor Capacity: 15mlx6 Max Speed: 4000r/min Max RCF: 1980×g</p>
Model	AE-TD4Y																						
Max Speed	4000r/min																						
Max RCF	1980×g																						
Max Capacity	15mlx6																						
Speed Accuracy	±30r/min																						
Time Setting	1min-99min																						
Noise	<65dB(A)																						
Power Supply	110V/60Hz, 220V/50Hz																						
Dimension(LxWxH)	310mm×260mm×220mm																						
Packing Size(LxWxH)	370mm×320mm×280mm																						
Net Weight	14kg																						



# TECHNOCRATIC SERVICES

ISO 9001:2015 Certified Company



**Please feel free to contact Ms Erum manzoor 021-35001798,  
0333-0297007, 0301-9228811.**

#### Head Office Karachi

- +92-21-35001798
- +92-333-0297007
- +92-321-2202241
- +92-300-3972448
- www.technocratics.com.pk
- technoratics@gmail.com
- info@technocratics.com.pk
- Plot No. 763, Future Colony Mor, Main Korangi Industrial Area Road, Karachi-75120

#### Branch Office Lahore

- + 92 301 9228811, + 92 334 6600557
- L-32 Main Street, Al-Nafy Homes, Ring Road, Near WAPDA Office Harbanspura, Lahore