




## Mini Centrifuge

Description	Actual picture														
<p><b>Mini Centrifuge / Micro Centrifuge</b></p> <p><b>Model # MC-4/6</b></p> <p>MC-4/6 Series Mini-Centrifuge is compact and cost effective, allowing any workstation to be equipped with a "personal" mini centrifuge. It includes 2 quick release interchangeable rotors, two speed options (4,000rpm &amp; 6000rpm). The electronic brake can reach quick deceleration to reduce handling time.</p> <p><b>Technical Specification</b></p> <table border="1" data-bbox="272 968 1081 1497"> <tr> <td><b>Model</b></td> <td>MC-4/6-3</td> </tr> <tr> <td><b>Capacity</b></td> <td>48well 50ul PCR plate rotor 4x5ml Eppendorf-tubes rotor</td> </tr> <tr> <td><b>Two Speed Options</b></td> <td>4,000rpm &amp; 6,000rpm</td> </tr> <tr> <td><b>Max. RCF</b></td> <td>1,000xg (4000rpm) 2,000xg (6000rpm)</td> </tr> <tr> <td><b>Dimensions</b></td> <td>135x147x107mm</td> </tr> <tr> <td><b>Weight</b></td> <td>1kg</td> </tr> <tr> <td><b>Power Supply</b></td> <td>AC 100-240V, 50/60Hz</td> </tr> </table> <p><b>Make :</b> A&amp;E LAB UK Co, Ltd</p>	<b>Model</b>	MC-4/6-3	<b>Capacity</b>	48well 50ul PCR plate rotor 4x5ml Eppendorf-tubes rotor	<b>Two Speed Options</b>	4,000rpm & 6,000rpm	<b>Max. RCF</b>	1,000xg (4000rpm) 2,000xg (6000rpm)	<b>Dimensions</b>	135x147x107mm	<b>Weight</b>	1kg	<b>Power Supply</b>	AC 100-240V, 50/60Hz	
<b>Model</b>	MC-4/6-3														
<b>Capacity</b>	48well 50ul PCR plate rotor 4x5ml Eppendorf-tubes rotor														
<b>Two Speed Options</b>	4,000rpm & 6,000rpm														
<b>Max. RCF</b>	1,000xg (4000rpm) 2,000xg (6000rpm)														
<b>Dimensions</b>	135x147x107mm														
<b>Weight</b>	1kg														
<b>Power Supply</b>	AC 100-240V, 50/60Hz														



# TECHNOCRATIC SERVICES

ISO 9001:2015 Certified Company



**Mini Centrifuge / Micro Centrifuge (Model # MC-4/6)**

#### Head Office Karachi

+92-21-35001798  
+92-333-0297007  
+92-321-2202241  
+92-300-3972448

[www.technocratics.com.pk](http://www.technocratics.com.pk)  
[technoratics@gmail.com](mailto:technoratics@gmail.com)  
[info@technocratics.com.pk](mailto:info@technocratics.com.pk)

Plot No. 763, Future Colony Mor, Main Korangi  
Industrial Area Road, Karachi-75120

#### Branch Office Lahore

+ 92 301 9228811, + 92 334 6600557

L-32 Main Street, Al-Nafy Homes, Ring Road,  
Near WAPDA Office Harbanspura, Lahore