








**TECHNOCRATIC SERVICES**

ISO 9001:2015 Certified Company

## Wipes for Clean room/ sterile are, Inject able Are for Class 100-10000

Description	Image
<p><b><u>Wipes for clean room/Sterile area /Inject able area(Lint Free):</u></b></p> <p><b>Features:</b></p> <ul style="list-style-type: none"> <li>- Laser sealed edge ensures ultra low lint and particle level</li> <li>- This package should therefore be opened in a clean room only</li> <li>- Good absorbency and wiping efficiency</li> <li>- Resists abrasion under rigorous use</li> <li>- Customizable size, packaging and grammage</li> <li>- wipe sensitive surface without any damage</li> </ul> <p><b>Specifications:</b></p> <p>Model Number: M-3</p> <p>Material: Non-woven Fabric, 100%(not lint free)</p> <p>Color: White</p> <p>Weight: 35gsm</p> <p>Size: 9" x 9"</p> <p>Style: Four folds</p> <p>Cutting way: Laser cut, cold cut, ultrasonic cut</p> <p>Packing: 100pcs/bag,</p> <p>Certificate: ISO</p> <p>Make: Jiangsu, China</p>	  
<p><b><u>Wipes for clean room/Sterile area /Inject able area(Lint Free):</u></b></p> <p><b>Features:</b></p> <ul style="list-style-type: none"> <li>- Laser sealed edge ensures ultra low lint and particle level</li> <li>- This package should therefore be opened in a clean room only</li> <li>- Good absorbency and wiping efficiency</li> <li>- Resists abrasion under rigorous use</li> <li>- Customizable size, packaging and grammage</li> <li>- wipe sensitive surface without any damage</li> </ul> <p><b>Specifications:</b></p> <p>Material: Nonwoven</p> <p>Brand Name: YCtek</p> <p>Material: Spunlaced Nonwoven, 55% wood pulp &amp; 45% dry polyester</p> <p>Weight: 50-75gsm</p> <p>Size: 9" x 9"</p> <p>Feature: Anti-Static, Eco-Friendly, Shrink-Resistant, Tear-Resistant</p> <p>Certification: ISO, RoHS, SGS</p> <p>Pattern: Plain</p> <p>Features: Low lint, high absorbency, softness</p> <p>Packing: 300 pcs/bag</p> <p>Color: White</p> <p>Make: Jiangsu, China</p>	 



# TECHNOCRATIC SERVICES

ISO 9001:2015 Certified Company

## Wipes for clean room/Sterile area /Inject able area(Lint Free):

### Features:

- Laser sealed edge ensures ultra low lint and particle level
- This package should therefore be opened in a clean room only
- Good absorbency and wiping efficiency
- Resists abrasion under rigorous use
- Customizable size, packaging and grammage
- wipe sensitive surface without any damage

### Specifications:

Item:	Wipes
Material:	Polyester
Color:	White
Product name:	Microfiber cleaning cloth clean room wiper
Application:	Laboratory, Sterile are, clean room
Feature:	Eco-Friendly, dust free, Lint free
Usage:	Clean room Protection
Type:	clean room wiper electronic wiper
Weight:	140gm
Cut Edge:	Laser cut or Ultrasonic cut
Clean Class:	100—10000
Packaging:	100 pcs/Bag
Make:	SHENZHEN, China



#### Head Office Karachi



#### Branch Office Lahore

