



Vertical autoclave high pressure steam sterilizer 35 liter

Description	Image
<p>Product description: MODEL: LS-35LD</p> <p>It is installed with an electrical heater, an automatically pressure-temperature controller, a safety valve, a releasing valve, a pressure-temperature indicator, an alarm bell for finished sterilization and instrument to automatically cut off the power of the heater. It has the advantages of effective sterilization, convenient operation, safe, less electro city consuming and cheap. It is ideal equipment for clinics, scientific research institutes and other organizations to sterilizer surgical instruments, fabric, glassed, culture media etc.</p> <p>Main features:</p> <ol style="list-style-type: none">1:Fully stainless high-alloy steel structure2:Digital display of working status, touch of key3:Manual or semi-automatic control4:Auto discharge the cool air, and steam discharging automatically after sterilization5:Automatically shut off with beep reminding after sterilization6:The accident protection system7:Automatic off of heating elements at reducing of a water level8:Blocking top cover at working9:Blockage of starting at opening10:Heater protection sensor11:Pressure and water level control12:With two stainless steel sterilizing baskets13:Easy to operate, safe and reliable14:The drying system is optional, which can be equipped according to the customers request15:Printer system could be install according to request.	  



TECHNOCRATIC SERVICES

ISO 9001:2015 Certified Company

Chamber volume	35 L
Working pressure	0.1 – 0.22 Mpa
Working temperature	120 – 134 C
Max working pressure	0.23 Mpa
Heat average	≤ 0.1C
Timer	0-99 min or 0 -99 hr 59 minute
Adjustment of temp	105 – 134 C
Power	2.5 kW/AC 220V 50Hz
Over all dimension	480*460*850mm
G.Wt/Net/wt	56kg/42kg



Head Office Karachi

- +92-21-35001798
- +92-333-0297007
- +92-321-2202241
- +92-300-3972448
- www.technocratics.com.pk
- technoratics@gmail.com
- info@technocratics.com.pk
- Plot No. 763, Future Colony Mor, Main Korangi Industrial Area Road, Karachi-75120

Branch Office Lahore

- 0333-0297007, 0335-4944616, 0321-2202241
- Plot # 42, Street # 2, Block-C, Asim Town, Harbanspura Underpass Near Ring Road, Lahore-Pakistan.