






Vertical autoclave high pressure steam sterilizer 50 liter

Description	Image																		
<p>Vertical Autoclave (Steam Sterilizer) Model: SSV50HD</p> <p>It is installed with an electrical heater, an automatically pressure-temperature controller, a safety valve, a releasing valve, a pressure-temperature indicator, an alarm bell for finished sterilization and instrument to automatically cut off the power of the heater. It has the advantages of effective sterilization, convenient operation, safe, less electricity consuming and cheap. It is ideal equipment for clinics, scientific research institutes and other organizations to sterilize surgical instruments, fabric, glassed, culture media etc.</p> <p>Main features:</p> <ol style="list-style-type: none"> 1: Fully stainless high-alloy steel structure 2: Digital display of working status, touch of key 3: Manual or semi-automatic control 4: Auto discharge the cool air, and steam discharging automatically after sterilization 5: Automatically shut off with beep reminding after sterilization 6: The accident protection system 7: Automatic off of heating elements at reducing of a water level 8: Blocking top cover at working. 9: Blockage of starting at opening 10: Heater protection sensor 11: Pressure and water level control 12: With two stainless steel sterilizing baskets 13: Easy to operate, safe and reliable. <p>Technical Specifications:</p> <table border="0"> <tr> <td>Instrument classification:</td> <td>Class II</td> </tr> <tr> <td>Chamber volume 50L:</td> <td>Working pressure: 0.22MPa</td> </tr> <tr> <td>Working temperature:</td> <td>120-134°C</td> </tr> <tr> <td>Max. working pressure:</td> <td>0.23 Mpa</td> </tr> <tr> <td>Adjustment of temperature:</td> <td>105-134°C</td> </tr> <tr> <td>Timer :</td> <td>0-60min</td> </tr> <tr> <td>LCD display of working status:</td> <td>touch type key.</td> </tr> <tr> <td>Hand wheel type open door structure:</td> <td>movable type</td> </tr> <tr> <td>Place of Origin:</td> <td>Shanghai, China</td> </tr> </table>	Instrument classification:	Class II	Chamber volume 50L:	Working pressure: 0.22MPa	Working temperature:	120-134°C	Max. working pressure:	0.23 Mpa	Adjustment of temperature:	105-134°C	Timer :	0-60min	LCD display of working status:	touch type key.	Hand wheel type open door structure:	movable type	Place of Origin:	Shanghai, China	  
Instrument classification:	Class II																		
Chamber volume 50L:	Working pressure: 0.22MPa																		
Working temperature:	120-134°C																		
Max. working pressure:	0.23 Mpa																		
Adjustment of temperature:	105-134°C																		
Timer :	0-60min																		
LCD display of working status:	touch type key.																		
Hand wheel type open door structure:	movable type																		
Place of Origin:	Shanghai, China																		

Head Office Karachi

-  +92-21-35001798
-  www.technocratics.com.pk
-  +92-333-0297007
-  technoratics@gmail.com
-  +92-321-2202241
-  info@technocratics.com.pk
-  Plot No. 763, Future Colony Mor, Main Korangi Industrial Area Road, Karachi-75120

Branch Office Lahore

-  0333-0297007, 0335-4944616, 0321-2202241
-  Plot # 42, Street # 2, Block-C, Asim Town, Harbanspura Underpass Near Ring Road, Lahore-Pakistan.