



# TECHNOCRATICS (Pvt) Ltd.

ISO 9001:2015 Certified Company  
NTN # A351529 GST # 3277876261645



## Vacuum Dust collector for capsule filling machine, capsule polish, tablet compression, bottle or blowing or carton de-dusting

### Vacuum Dust collector (For capsule filling , tablet compression, bottle or blowing or carton de-dusting)

#### Model XBG-1100

- Large exhaust air volume, strong vacuum suction;
- High quality stainless steel structure;
- Filter bag adopts manual vibration dust removal;
- Simple operation and easy maintenance;
- Smooth operation, low noise.

#### Application

Stainless steel SS304 vacuum chamber suitable for all kinds of high vacuum and super-high vacuum condition. It is (For capsule filling , tablet compression, bottle or blowing or carton de-dusting)

And also widely used in optical coating, PVD, aerospace industry, semi-conductor industry etc.

Powder can be recovered and it conforms to GMP standard. Silencer-device, deduce the noise.

Rated voltage	380V
Power	1100W
Maximum flow rate (m/h)	135
Noise (DB)	≤60
Maximum voltage (kpa)	≤12
Vacume degree (kpa)	15
Overall dimensions(mm)	700x450x1200
Vibration	Manual
Filtration	Single bag
Make:	Pharmaland, Taiwan

Rated voltage	380V
Power	2200W
Maximum flow rate (m/h)	210
Noise (DB)	≤60
Maximum voltage (kpa)	≤18
Vacume degree (kpa)	20
Overall dimensions(mm)	850x500x1300
Make:	Pharmaland, Taiwan



**Head Office Karachi**  
+92-21-35001798  
+92-333-0297007  
+92-321-2202241  
+92-300-3972448  
www.technocratics.com.pk  
technocratics@gmail.com  
info@technocratics.com.pk  
Plot No. 763, Future Colony Mor, Main Korangi Industrial Area Road, Karachi-75120

**Branch Office Lahore**  
+92-301-9228811, +92-334-6600557  
L-32 Main Street, Al-Nafy Homes, Ring Road, Near WAPDA Office, Harbanspura, LAHORE  
0333-0297007, 0321-2202241

**Branch Office Islamabad**  
+92 334 6600557  
Awan Plaza River Garden, Islamabad Express Highway, Islamabad  
0333-0297007, 0321-2202241



# TECHNOCRATICS (Pvt) Ltd.




ISO 9001:2015 Certified Company  
NTN # A351529 GST # 3277876261645



## Head Office Karachi

+92-21-35001798  
+92-333-0297007  
+92-321-2202241  
+92-300-3972448  
 [www.technocratics.com.pk](http://www.technocratics.com.pk)  
 [technocratics@gmail.com](mailto:technocratics@gmail.com)  
[info@technocratics.com.pk](mailto:info@technocratics.com.pk)  
Plot No. 763, Future Colony Mor, Main Korangi  
Industrial Area Road, Karachi-75120

## Branch Office Lahore

 +92-301-9228811, +92-334-6600557  
 L-32 Main Street, Al-Nafy Homes, Ring Road,  
Near WAPDA Office, Harbanspura, LAHORE  
 0333-0297007, 0321-2202241

## Branch Office Islamabad

 +92 334 6600557  
 Awan Plaza River Garden,  
Islamabad Express Highway, Islamabad  
 0333-0297007, 0321-2202241