



TECHNOCRATICS (Pvt) Ltd.

ISO 9001:2015 Certified Company

Wipes for clean room/Sterile area /Inject able area(Lint Free):

Features:

- Laser sealed edge ensures ultra low lint and particle level.
- This package should therefore be opened in a clean room only
- Good absorbency and wiping efficiency.
- Resists abrasion under rigorous use.
- Customizable size, packaging and grammage.
- wipe sensitive surface without any damage.

Specifications:

Sr. #	Technical Parameter:	
1	Material	Nonwoven
2	Brand Name	YCtek
3	Material	Spun laced Nonwoven, 55% wood pulp & 45% dry polyester
4	Weight	50-75gsm
5	Size	9" x 9"
6	Feature	Anti-Static, Eco-Friendly, Shrink-Resistant, Tear-Resistant
7	Certification	ISO, RoHs, SGS
8	Pattern	Plain
9	Features	Low lint, high absorbency, softness
10	Packing	300 pcs/bag
11	Color	White
10	Make	Jiangsu, China

Picture:



Head Office Karachi

+92-21-35001798
+92-333-0297007
+92-321-2202241
+92-300-3972448
www.technocratics.com.pk
technocratics@gmail.com
info@technocratics.com.pk
Plot No. 763, Future Colony Mor, Main Korangi Industrial Area Road, Karachi-75120

Branch Office Lahore

+92-301-9228811, +92-334-6600557
L-32 Main Street, Al-Nafy Homes, Ring Road, Near WAPDA Office, Harbanspura, LAHORE
0333-0297007, 0321-2202241

Branch Office Islamabad

+92 334 6600557
Awan Plaza River Garden, Islamabad Express Highway, Islamabad
0333-0297007, 0321-2202241



TECHNOCRATICS (Pvt) Ltd.

ISO 9001:2015 Certified Company

--	--

Head Office Karachi

+92-21-35001798
+92-333-0297007
+92-321-2202241
+92-300-3972448
www.technocratics.com.pk
technocratics@gmail.com
info@technocratics.com.pk
Plot No. 763, Future Colony Mor, Main Korangi
Industrial Area Road, Karachi-75120

Branch Office Lahore

+92-301-9228811, +92-334-6600557
L-32 Main Street, Al-Nafy Homes, Ring Road,
Near WAPDA Office, Harbanspura, LAHORE
0333-0297007, 0321-2202241

Branch Office Islamabad

+92 334 6600557
Awan Plaza River Garden,
Islamabad Express Highway, Islamabad
0333-0297007, 0321-2202241