



TECHNOCRATICS (Pvt) Ltd.

ISO 9001:2015 Certified Company
NTN # A351529 GST # 3277876261645



Vacuum Dust collector for capsule filling machine, capsule polish, tablet compression, bottle or blowing or carton de-dusting

Vacuum Dust collector

(For capsule filling , tablet compression, bottle or blowing or carton de-dusting)

Model XBG-1100/2200W

- Large exhaust air volume, strong vacuum suction;
- High quality stainless steel structure;
- Filter bag adopts manual vibration dust removal;
- Simple operation and easy maintenance;
- Smooth operation, low noise.

Application

Stainless steel SS304 vacuum chamber suitable for all kinds of high vacuum and super-high vacuum condition. It is (For capsule filling , tablet compression, bottle or blowing or carton de-dusting)

And also widely used in optical coating, PVD, aerospace industry, semi-conductor industry etc.

Powder can be recovered and it conforms to GMP standard. Silencer-device, deduce the noise.

Rated voltage	380V
Power	1100W
Maximum flow rate (m/h)	135
Noise (DB)	≤60
Maximum voltage (kpa)	≤12
Vacume degree (kpa)	15
Overall dimensions(mm)	700x450x1200
Vibration	Manual
Filtration	Single bag
Make:	Pharmaland, Taiwan

Rated voltage	380V
Power	2200W
Maximum flow rate (m/h)	210
Noise (DB)	≤60
Maximum voltage (kpa)	≤18
Vacume degree (kpa)	20
Overall dimensions(mm)	850x500x1300
Make:	Pharmaland, Taiwan



Head Office Karachi
+92-21-35001798
+92-333-0297007
+92-321-2202241
+92-300-3972448
www.technocratics.com.pk
technocratics@gmail.com
info@technocratics.com.pk
Plot No. 763, Future Colony Mor, Main Korangi Industrial Area Road, Karachi-75120

Branch Office Lahore
+92-301-9228811, +92-334-6600557
L-32 Main Street, Al-Nafy Homes, Ring Road, Near WAPDA Office, Harbanspura, LAHORE
0333-0297007, 0321-2202241

Branch Office Islamabad
+92 334 6600557
Awan Plaza River Garden, Islamabad Express Highway, Islamabad
0333-0297007, 0321-2202241



TECHNOCRATICS (Pvt) Ltd.

ISO 9001:2015 Certified Company
NTN # A351529 GST # 3277876261645



Head Office Karachi

+92-21-35001798
+92-333-0297007
+92-321-2202241
+92-300-3972448
 www.technocratics.com.pk
 technocratics@gmail.com
info@technocratics.com.pk
Plot No. 763, Future Colony Mor, Main Korangi
Industrial Area Road, Karachi-75120

Branch Office Lahore

 +92-301-9228811, +92-334-6600557
 L-32 Main Street, Al-Nafy Homes, Ring Road,
Near WAPDA Office, Harbanspura, LAHORE
 0333-0297007, 0321-2202241

Branch Office Islamabad

 +92 334 6600557
 Awan Plaza River Garden,
Islamabad Express Highway, Islamabad
 0333-0297007, 0321-2202241