



TECHNOCRATICS (Pvt) Ltd.

ISO 9001:2015 Certified Company
NTN # A351529 GST # 3277876261645



Sodium Lamp for polarimeter

S. No	Description	Image								
01	Sodium Lamp Special purpose sodium lamp for WXG-4 manual polarimeter. Description: • Monochromatic sodium light source • Wavelength: 589.4 nm • Power requirement: 220V/50Hz Technical specifications <table><tr><td>Material</td><td>Glass and metal</td></tr><tr><td>Connector</td><td>4-pin</td></tr><tr><td>Dimensions</td><td>155 × 28(Ø) mm</td></tr><tr><td>Weight</td><td>70 g</td></tr></table>	Material	Glass and metal	Connector	4-pin	Dimensions	155 × 28(Ø) mm	Weight	70 g	
Material	Glass and metal									
Connector	4-pin									
Dimensions	155 × 28(Ø) mm									
Weight	70 g									

For: Technocratics , Erum Manzoor 03330297007 ,021-35001798 ,0300-3972448