







TECHNOCRATICS (Pvt) Ltd.

ISO 9001:2015 Certified Company
NTN # A351529 GST # 3277876261645



Suture Needle Tube Toughness Testing Machine

Description		Image
Suture Needle Tube Toughness Testing Machine toughness measuring instrument. Suitable for infusion set, blood transfusion set, infusion needle, anesthesia kit, pipeline, catheter, filter, quick connector, etc. According to different uses and designs, infusion needles can be divided into intravenous needles, arterial needles, local injection needles and so on.		
Injection Needle Tip Needle Toughness Tester/Suture Needle Tube Toughness Testing Machine		
Name	Injection needle toughness tester	
Test frequency	0.5 HZ	
Bending Angle	15°, 20°, 25(25°4590°)3 positions, switch set bending Angle, 4 high brightness LED display;	
Bending distance	855mm	
Accuracy	±0.1 mm	
Number of experiments	Bidirectional force applied 20 times and fatigue fracture test	
Number of experiments	Bidirectional force applied 20 times and fatigue fracture test	
Net weight	10kg	
Overall dimensions	Length * Width * Height :500X560X540(mm)	
Scope of application	It is specially used to test the toughness of medical indwelling n mechanism swing mechanism, fixture and other parts. The intelligent controller consists of CPU, motherboard, LCD display, power supply and other parts	

Contact: 03330297007
ERUM MAZOOR



TECHNOCRATICS (Pvt) Ltd.

ISO 9001:2015 Certified Company
NTN # A351529 GST # 3277876261645



Head Office Karachi



+92-21-35001798



www.technocratics.com.pk



+92-333-0297007



technocratics@gmail.com

+92-321-2202241

+92-300-3972448

info@technocratics.com.pk



Plot No. 763, Future Colony Mor, Main Korangi
Industrial Area Road, Karachi-75120

Branch Office Lahore



0333-0297007, 0335-4944616, 0321-2202241



Plot # 42, Street # 2, Block-C, Asim Town,
Harbanspura Underpass Near Ring Road,
Lahore-Pakistan.