



TECHNOCRATICS (Pvt) Ltd.

ISO 9001:2015 Certified Company

R & D Tablet Compression Machine

SINGLE PUNCH TABLET PRESS MACHINE WITH 01 ROUND DIE FREE

Description

SINGLE PUNCH TABLET PRESS MACHINE :

With 01 round die free

Model No: TDP-0

- DP 0 Tablet Press is an advanced, manually operated pill press that is able to press round tablets up to 10 mm in diameter and 6 mm thick without the need for electrical power.
- Small and light weight, making it far more portable than electric powered tablet presses.
- Easy to use and maintain
- Interchangeable dies, No running costs
- Contact parts Stainless steel

Technical Specification:

Number of dies	1
Production capacity	30-50/min
Diameter of the tablets	5-10 mm
Thickness of tablet	2-6 mm
Max pressure	3.1 kN
Number of filling stations	1
Double layered tablet	No
Weight	25.4 kg / 56 lbs
Height	51 cm / 20 in
Width	30 cm / 11.8 in
Depth	22 cm / 8.6 in

Make. Hong Kong, China

Head Office Karachi

+92-21-35001798
+92-333-0297007
+92-321-2202241
+92-300-3972448
www.technocratics.com.pk
technocratics@gmail.com
info@technocratics.com.pk
Plot No. 763, Future Colony Mor, Main Korangi Industrial Area Road, Karachi-75120

Branch Office Lahore

+92-301-9228811, +92-334-6600557
L-32 Main Street, Al-Nafy Homes, Ring Road, Near WAPDA Office, Harbanspura, LAHORE
0333-0297007, 0321-2202241

Branch Office Islamabad

+92 334 6600557
Awan Plaza River Garden, Islamabad Express Highway, Islamabad
0333-0297007, 0321-2202241



TECHNOCRATICS (Pvt) Ltd.

ISO 9001:2015 Certified Company



Head Office Karachi

+92-21-35001798
+92-333-0297007
+92-321-2202241
+92-300-3972448
www.technocratics.com.pk
technocratics@gmail.com
info@technocratics.com.pk
Plot No. 763, Future Colony Mor, Main Korangi Industrial Area Road, Karachi-75120

Branch Office Lahore

+92-301-9228811, +92-334-6600557
L-32 Main Street, Al-Nafy Homes, Ring Road,
Near WAPDA Office, Harbanspura, LAHORE
0333-0297007, 0321-2202241

Branch Office Islamabad

+92 334 6600557
Awan Plaza River Garden,
Islamabad Express Highway, Islamabad
0333-0297007, 0321-2202241