



TECHNOCRATICS (Pvt) Ltd.

ISO 9001:2015 Certified Company
NTN # A351529 GST # 3277876261645



Sodium lamp for polarimeter

S. No	Description	Image																		
01	Polarimeter Tubes Glass polarimeter tubes for manual, semiautomatic and automatic polarimeters. Description: • General purpose polarimeter tube with bubble trap • Glass body, metal collars • Selectable lengths are 100 and 200 mm Technical specifications: <table> <tr> <td>Material</td><td>Glass and metal</td><td>Glass and metal</td></tr> <tr> <td>Glass Body Length</td><td>100 mm</td><td>200 mm</td></tr> <tr> <td>Glass Body Diameter</td><td>15 mm</td><td>15 mm</td></tr> <tr> <td>Sample Volume</td><td>Approximately 20 mL</td><td>Approximately 40 mL</td></tr> <tr> <td>External Diameter</td><td>27 mm</td><td>27 mm</td></tr> <tr> <td>Weight</td><td>200 g</td><td>250 g</td></tr> </table>	Material	Glass and metal	Glass and metal	Glass Body Length	100 mm	200 mm	Glass Body Diameter	15 mm	15 mm	Sample Volume	Approximately 20 mL	Approximately 40 mL	External Diameter	27 mm	27 mm	Weight	200 g	250 g	 
Material	Glass and metal	Glass and metal																		
Glass Body Length	100 mm	200 mm																		
Glass Body Diameter	15 mm	15 mm																		
Sample Volume	Approximately 20 mL	Approximately 40 mL																		
External Diameter	27 mm	27 mm																		
Weight	200 g	250 g																		

For: Technocratics , Erum Manzoor 03330297007 ,021-35001798 ,0300-3972448