

NITRILE RUBBER LONG CUFF GLOVES

Picture Description:

The image features the Nitrile Rubber Long Cuff Gloves, showing the gloves designed for work and safety sterilization, highlighting their long cuff for extended forearm protection.

Technical Data:

Product Name: Nitrile Rubber Long Cuff Gloves

Material: Nitrile Rubber (Latex-Free)

Color: White / Off-White

Weight: 190 g

Length: Approx. 80 cm (Long cuff design)

Cuff Width: 24 ± 2 cm

Thickness: 0.9 ± 0.1 mm

Surface: Smooth / Slightly Textured

Sterility: Non-sterile / Sterilizable

Powder: Powder-Free

Sizes Available: S, M, L, XL

Resistance: Chemical, Oil, Puncture Resistant

Cuff Type: Extended Forearm Protection

Packaging: Sealed Bag

Applications: Glove box operations, sterilization testing, pharmaceutical testing, chemical industry, etc.

Key Features:

- Latex-Free Material: Safe for users with latex allergies.
- Extra-Long Cuff Design: Provides full forearm protection (approx. 80 cm).
- Excellent Chemical Resistance: Suitable for harsh environments.
- High Durability: Resistant to tearing and punctures.
- Comfortable Fit: Flexible and ergonomic for long use.
- Powder-Free: Reduces contamination risk.
- Sterilization Compatible: Can be used in controlled environments.
- Good Grip Control: Ensures safe handling of tools and materials.

Applications:

- Pharmaceutical & Laboratory Work
- Food Processing & Handling
- Chemical Handling & Industrial Use
- Cleaning & Sanitation Operations
- Medical & Healthcare (Non-Surgical)
- Agriculture & Veterinary Work
- Automotive & Industrial Work
- Sterilization & Decontamination Tasks

NITRILE RUBBER LONG CUFF GLOVES

Size Chart:

- **S:** Palm Width: 7.5 – 8.0 cm, Cuff Width: 24 ± 2 cm
- **M:** Palm Width: 8.5 – 9.0 cm, Cuff Width: 24 ± 2 cm
- **L:** Palm Width: 9.0 – 9.5 cm, Cuff Width: 24 ± 2 cm
- **XL:** Palm Width: 9.5 – 10.0 cm, Cuff Width: 24 ± 2 cm

